Accutrend Control CH 1, Roche Diagnostics

Accutrend Control CH 1, Roche Diagnostics

Accutrend® Control CH 1

Accutrend Control CH 1, Roche Diagnostics

Accutrend Control CH 1, Roche Diagnostics

DE

IN VITRO DIAGNOSTICUM

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ES

FOR IN VITRO DIAGNOSTIC USE

For performance checks of the Accutrend GCT or Accutrend Plus meter using Accutrend Cholesterol test strips. Suitable for self-testing.

General information : Use Accutrend Control CH 1 control solution to check the performance of your Accutrend GC, Accutrend GCT or Accutrend Plus meter and the Accutrend Cholesterol test strips and their correct handling. Perform a performance check in the following situations:

- When learning how to measure cholesterol
- after changing batteries
- after opening a new pack of test strips and re-calibrating your meter
- after cleaning the meter
- if you have reasons to doubt the validity of a reading

How the test works: Each test strip has a test area containing detection reagents. When control solution is applied, a chemical reaction takes place and the test area changes colour. The meter records this change in colour and converts the measurement signal to the displayed result using the data previously entered via the code strip.

Important notes: Do not perform performance checks at temperatures below +18°C or above +30°C. Allow test strips and control solution just removed from the refrigerator to reach a temperature within this range before use. Use a fresh test strip if the drop of solution has not completely covered the test area, if several drops were applied or if a drop contained air bubbles. Do not touch the tip of the bottle.

Storage and stability: The expiry date printed on the box and the label is valid when stored between +2°C and +25°C. Do not freeze. When stored at the conditions stated above, opened bottles of cholesterol control solution are stable for 3 months but not beyond the the expiry date given on the pack. When opening the bottle for the first time, write the date onto the bottle label. Disposal: Empty the bottle and discard as regular packaging waste.

Sources of error: If a performance check with the Accutrend Control CH 1 control solution produces a reading outside the stated range, please check:

- that the cholesterol measurement was carried out correctly (see operating instructions)
- that the test strip guide is clean (follow instructions for use)
- that the test strips and the control solution used were correctly stored (expiry date, storage conditions)
- that meter, test strips and control solution were stored at the correct temperature range of +18°C bis +30°C during the measurement.

Preparing and performing the performance check: When carrying out a performance check, please follow the instructions given in this insert, in the operating instructions of your meter and in the package insert of the test strip. When carrying out a performance check, proceed in the same way as you do for a cholesterol measurement; however, instead of blood use Accutrend Control CH 1 control solution. Apply a drop of control solution to the test pad; make sure it is completely covered. Hold the dropper vertically with its tip showing down; do not touch the test pad with the tip. Wipe the tip with a clean paper towel or a swab before and after use and then close it again tightly. The acceptable control

EN

range for a given test strip lot (code number) in combination with the control solution lot is shown on one of the notes supplied with the test strip pack and the control solution pack. If the reading is within this range, both the meter and the test are working properly. If the reading is outside the range, check if all your handling steps were correct (see also „sources of error“). Eliminate all possible sources of error and repeat the performance check. If, in a second measurement, the reading is still outside the specified range, please contact your local Roche Diagnostics representative.

Contents of the pack: 1 bottle with cholesterol control solution, note.

Additionally required material: One of the above listed meters with operating instructions, Accutrend Cholesterol test strips, 1 paper towel.

Notes for doctors and laboratories:

Applications: Perform quality control steps according to guidelines established for your practice or laboratory.

Contents: Aqueous solution including chemical process materials and extract of animal origin. The origin of the biological component is as follows: protein/lipoprotein serum (bovine). The concentration of the components is lot-specific. The range of the component to be determined is listed in the note.

Performance characteristics of the Accutrend Control CH 1 control solution: Reproducibility (day-to-day imprecision): In typical measurement series, the majority of the measured coefficients of variation (CV) was < 3.8%.

Last revision: 2007-06

Distributed by: Roche Diagnostics Limited, Charles Avenue, Burgess Hill, West Sussex, RH15 9RY, **United Kingdom**

Roche Diagnostics Australia Pty Ltd., 31 Victoria Avenue, Castle Hill, NSW 2154, **Australia**

Roche Diagnostics, 201, Boulevard Armand-Frappier, Laval, Québec (CANADA) HZV 4A2, **Canada**

Roche Diagnostics N.Z. Ltd, 15 Rakino Way, Box 62-089, Mt. Wellington, Auckland, **New Zealand**

Presentation: Pack with 1.5 mL cholesterol control solution; good for approx. 30 performance checks. For technical reasons, the bottle is not filled to capacity. Accutrend Control CH 1, REF 1 1418289190.

For an explanation of the symbols used please refer to the end of this package insert.

EN

FR

DE

IT

ES

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE




DE

DE

DE

DE

DE

	Bitte Gebrauchsanweisung beachten / Lire attentivement les instructions d'utilisation / Consultare le istruzioni per l'uso / De gebruiksaanwijzing raadplegen, a.u.b. / Consulte las instrucciones de uso / Consultar as instruções de utilização / Please consult instructions for use
	Lagerung bei / Conserver entre / Conservare a / Bewaren bij / Conservare a / Conservar a / Store at
	Verwendbar bis / Date de péremption / Scadenza / Houdbaar to / Fecha de caducidad / Data limite de utilização / Use by / Expiry date
	Hersteller / Fabricant / Fabricante / Fabrikant / Fabbricante / Fabricante / Manufacturer
REF	Artikelnummer / Référence article / Codice dell'articolo / Número de Catálogo / Número de Referência / Catalogue number
LOT	Chargenbezeichnung / Numéro de lot / Numero di lotto / Lotnummer / Número del lote / Número do lote / Lot number
IVD	In-vitro-Diagnostikum / IN VITRO DIAGNOSTIC / Per uso diagnostico in vitro / Voor in vitro diagnostisch gebruik / Diagnóstico in vitro / Para utilização em diagnóstico in vitro / For in vitro diagnostic use
 0123	Dieses Produkt entspricht den Anforderungen der Richtlinie über In-vitro-Diagnostika 98/79/EG. / Ce produit répond aux exigences de la directive 98/79 CE relative aux dispositifs médicaux de diagnostic in vitro. / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro. / Dit product voldoet aan de richtlijn 98/79/EG inzake medische hulpmiddelen voor in-vitro diagnostiek. / Este producto cumple las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro. / Este produto cumpre os requisitos da Directiva Europeia n.º 98/79/CE para dispositivos médicos para diagnóstico in vitro. / This product fulfils the requirements of Directive 98/79/EC on in vitro diagnostic medical devices.
 0123	

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

DE

0 5053030001 (01)

R1 (black)