

Intended use

The Propulse® Ng Electronic Ear Irrigation System is a Class IIa medical device. Its intended use is for the removal of cerumen, keratin, or any other (non-hygroscopic) foreign body from the human ear by means of irrigation with water. The Propulse® Ng Electronic Ear Irrigation System is CE Marked and complies with

- Medical Devices Directive – 93/42/EEC
- Electrical Safety – EN60601-1-2
- Electro Magnetic Compatibility

Each Propulse® Ng ear irrigation system comprises:



- Propulse® Ng Unit
- Integrated QrX™ Handle
- Propulse® Footswitch
- 10 - Propulse® QrX™ Single-Use Tips
- 10 - Cleaning Tablets (via registration card)
- Power Supply
- Internal Rechargeable Power Pack
- User Manual

Propulse® Chlor Clean cleaning tablets

Cleaning tablets for your Propulse® Ear Irrigator – the importance of using the correct strength of cleaning solution cannot be overstated. A solution that is too strong will in time damage the Propulse®, and one that is too weak will fail to provide the correct level of cleaning and decontamination. Propulse® cleaning tablets have been formulated to provide the exact strength of cleaning solution that is required to clean your Propulse® effectively and safely.

Please note: Damage caused to your Propulse® Ng Ear Irrigator as a result of not using the recommended cleaning agent will invalidate your warranty.

