



Why use ear tips for single use only?

Sanibel Supply® recommends using disposable ear tips for single use only to stop ear infections from spreading

At Sanibel Supply® we think that most audiologists would agree with us that it is crucial to avoid inadvertent transmission of disease within the audiology clinic. The good news is that it is quite easy to ensure that ear infections are not spread from one patient to another; use a new ear tip for each patient, and you will significantly reduce or eliminate the incidences of cross contamination.

Cerumen may be an infectious agent

Cerumen is not considered an infectious agent unless contaminated with blood or mucus, and due to its color and viscosity, visual detection of blood or ear discharge contaminants may be difficult. Therefore, you should treat cerumen as a potentially infectious agent. It is especially important to be mindful of using new tips when testing children and those with weaker immune systems, such as the frail and elderly as they are more susceptible to infections.

Are disposable ear tips expensive?

Sanibel Supply® offers volume discounts on ear tips. This will often offset the time and cost involved in the ear tip cleaning process.

Be more cost effective in your clinic

By using a new ear tip pr. every patient you reach a seal much more reliable than with re-used ear tips. How often have you not been fumbling for minutes before you get a seal with a patient? Using the right ear tip and a fresh one pr. patient ensures you a reliable fit with every patient, meaning you can calculate less minutes for testing a patient and in most cases allow you to have an extra patient pr. day without having stressed out personnel. The cost of using one ear tip pr. patient is far less than what you earn on your next patient, plus you can promote your clinic as taking patient safety and hygiene seriously.

6 strong reasons why you should use a new ear tip after each patient:

1

You reduce the risk of cross contamination between your patients

Unfortunately, it is not possible to be 100% certain that complete decontamination has taken place when cleaning ear tips.

2

You minimize the risk of non-sealing ear tips

Ear tips degrade over time, which can make it more difficult to attain a good seal, hence leading to less accurate results.

3

Your patients avoid allergic reactions from cleaning agent residue

Residue from cleaning agents may cause allergic reactions in some patients.

4

Your patients will appreciate it

Patients are more comfortable knowing that new ear tips are being used.

5

It's so much easier

Cleaning ear tips takes time, so why not give your staff the time to do what really matters?




























6

It's so much faster

A fresh ear tip creates a seal quicker and more reliable, so you do not waste time performing a standard procedure.



Selecting the correct ear tip

Age	Ear tips inserted against ear canal	Ear tips inserted against entrance of ear canal			Ear tips for Echo-Screen®	Ear tips for AccuScreen®
0.5 y/o		 Green pediatric			 4.0 mm, ANA	 3.7 mm, AOT
					 4.5 mm, ANA	 4.0 mm, AOT
		 3-5 mm, flanged	 3-6 mm, green	 6 mm, green	 5.0 mm, ANA	 4.5 mm, AOT
	 7.5 mm, umbrella	 4-7 mm, flanged	 6 mm, yellow		 infant tree, ANA	 5.0 mm, AOT
	 9 mm, umbrella	 5-8 mm, flanged			 mushroom, ANA	 4-7 mm, AOT
	 10 mm, umbrella	 7 mm, mushroom	 7 mm, red	 7 mm, yellow		
	 11 mm, umbrella	 8 mm, mushroom	 8 mm, blue	 8 mm, red		
	 13 mm, umbrella	 9 mm, mushroom	 9 mm, green	 9 mm, blue		
		 10 mm, mushroom	 10 mm, yellow	 10 mm, green		
	 15 mm, umbrella	 11 mm, mushroom	 11 mm, red	 11 mm, yellow		
10 y/o						
Adults		 12 mm, mushroom	 12 mm, blue	 12 mm, red		
	 19 mm, umbrella	 13 mm, mushroom	 13 mm, green	 13 mm, blue		
		 14 mm, mushroom	 14 mm, yellow	 14 mm, green		
		 15 mm, mushroom	 16 mm, red	 16 mm, yellow		
	 22 mm, umbrella	 19 mm, mushroom	 18.5 mm, blue			

Selecting the correct size of ear tip takes practice

In general, flanged or mushroom-shaped ear tips will be used for most testing (Impedance, DPOAE, TEOAE and Automated ABR).

Sometimes, an umbrella-shaped ear tip will be used for performing a quick impedance screening.

The table can be used as a guide to help you select the appropriate ear tip size.

Disposable Sanibel™ ADI ear tips

For use with the following acoustic immittance instruments and otoacoustic emission systems with CLEAR probe tip: **Amplivox**: all Otowave versions. **Interacoustics**: Titan, AT235, AA220, Eclipse, MT10, Sera, OtoRead II. **Maico**: easyTymp, touchTymp MI 22, MI 23, MI 24, MI 26, MI 34, MI 44, easyScreen, EroScan Plus,



ADI Assortment Box

Part no.	Description
8511445	Various sizes for Tympanometry and OAE
8501572	Screening ear tip kit, various sizes for newborn screening



Flanged infant ear tips

Part no.	Description	Unit
8502837	Pediatric (green)	Bag/25 pcs
8502838	Pediatric (green)	Bag/100 pcs
8012963	3-5 mm (red)	Bag/25 pcs
8012964	3-5 mm (red)	Bag/100 pcs
8012965	4-7 mm (blue)	Bag/25 pcs
8012966	4-7 mm (blue)	Bag/100 pcs
8012967	5-8 mm (yellow)	Bag/25 pcs
8012968	5-8 mm (yellow)	Bag/100 pcs

VOLUME DISCOUNT AVAILABLE for minimum orders of 10 bags



VOLUME DISCOUNT AVAILABLE for minimum orders of 10 bags



Medical grade silicone rubber
Biocompatible



umbrella ear tips

Part no.	Description	Unit
8012985	7.5 mm (green)	Bag/25 pcs
8012986	7.5 mm (green)	Bag/100 pcs
8012987	9 mm (red)	Bag/25 pcs
8012988	9 mm (red)	Bag/100 pcs
8012989	10 mm (blue)	Bag/25 pcs
8012990	10 mm (blue)	Bag/100 pcs
8012991	11 mm (yellow)	Bag/25 pcs
8012992	11 mm (yellow)	Bag/100 pcs
8012993	13 mm (green)	Bag/25 pcs
8012994	13 mm (green)	Bag/100 pcs
8012995	15 mm (red)	Bag/25 pcs
8012996	15 mm (red)	Bag/100 pcs
8012997	19 mm (blue)	Bag/25 pcs
8012998	19 mm (blue)	Bag/100 pcs
8012999	22 mm (yellow)	Bag/25 pcs
8013000	22 mm (yellow)	Bag/100 pcs



VOLUME DISCOUNT AVAILABLE for minimum orders of 10 bags



Mushroom ear tips

Part no.	Description	Unit
8013001	7 mm (blue)	Bag/25 pcs
8013002	7 mm (blue)	Bag/100 pcs
8013003	8 mm (yellow)	Bag/25 pcs
8013004	8 mm (yellow)	Bag/100 pcs
8012969	9 mm (green)	Bag/25 pcs
8012970	9 mm (green)	Bag/100 pcs
8012971	10 mm (red)	Bag/25 pcs
8012972	10 mm (red)	Bag/100 pcs
8012973	11 mm (blue)	Bag/25 pcs
8012974	11 mm (blue)	Bag/100 pcs
8012975	12 mm (yellow)	Bag/25 pcs
8012976	12 mm (yellow)	Bag/100 pcs
8012977	13 mm (green)	Bag/25 pcs
8012978	13 mm (green)	Bag/100 pcs
8012979	14 mm (red)	Bag/25 pcs
8012980	14 mm (red)	Bag/100 pcs
8012981	15 mm (blue)	Bag/25 pcs
8012982	15 mm (blue)	Bag/100 pcs
8012983	19 mm (yellow)	Bag/25 pcs
8012984	19 mm (yellow)	Bag/100 pcs



Single use ear tips for immittance instruments & OAE screeners

For use with many acoustic immittance instruments and otoacoustic emission systems, including the **MAICO** EroScan, **Interacoustics** OtoRead and **Grason-Stadler** TympStar.



VOLUME DISCOUNT AVAILABLE for minimum orders of 10 bags

Infant ear tips

Part no.	Description	Unit
8013005	3 mm (yellow)	Bag/25 pcs
8013006	3 mm (yellow)	Bag/100 pcs
8002036	4 mm (red)	Bag/25 pcs
8002035	4 mm flanged (red)	Bag/100 pcs



IA series ear tips

For use with many acoustic immittance instruments and otoacoustic emission systems, including: **GSI**: TymStar Pro, Allegro, Novus and Corti. **Interacoustics** (clear probe): AT235, AA222 and Titan and Otoread II, MT10. **Maico**: touchTym Racecar, touchTym MI24, MI26, MI34, MI36, easyTym. **Madsen/Otometrics**: Zodiac (new), Accuscreen. **Path Medical**: all Sentiero products.



Mushroom ear tips

Part no.	Description	Unit
8515491	Sample Kit (3-18.5mm), 3, 6-10 Count Each	
8520348	3-6 mm (green)	Bag/100 pcs
8520349	6 mm (yellow)	Bag/100 pcs
8520350	7 mm (red)	Bag/100 pcs
8520351	8 mm (blue)	Bag/100 pcs
8520352	9 mm (green)	Bag/100 pcs
8520353	10 mm (yellow)	Bag/100 pcs
8520354	11 mm (red)	Bag/100 pcs
8520355	12 mm (blue)	Bag/100 pcs
8520356	13 mm (green)	Bag/100 pcs
8520357	14 mm (yellow)	Bag/100 pcs
8520358	16 mm (red)	Bag/100 pcs
8520359	18.5 mm (blue)	Bag/100 pcs



ABR

IMP

OAE

IMPEDANCE EAR TIPS – Fading product

For use with the following systems with BLACK probe tip:
Interacoustics: Old MT10, AT235, AA22x, MI 22, MI 23, MI 24, MI 26, MI 34, MI 44. **FADING OUT** product: consider AKR instead.



VOLUME DISCOUNT AVAILABLE
for minimum orders of 10 bags

Umbrella ear tips

Part no.	Description	Unit
8107617	10 mm (green)	bag/100 pcs
8107618	13 mm (green)	bag/100 pcs
8107619	16 mm (green)	bag/100 pcs
8107620	19 mm (blue)	bag/100 pcs
8107621	22 mm (blue)	bag/100 pcs



IMP

AKR Mushroom ear tips

Use with the following acoustic immittance instruments and otoacoustic emission systems with BLACK probe tip from **Interacoustics**: AT235, AA220, AA220, Eclipse, MT10, Otoread (I), **Maico**: touchTym RaceCar MI 22, MI 23, MI 24, MI 26, MI 34, MI 44. **Amplivox**: Otowave (all with white probe). **Madsen/Otometrics**: Capella 2, Zodiac 901, Biologic AuDX Pro (+Flex). **Neurosoft**: Neuro-Audio. **Oscilla**: TSM 300/400 + T800 series. All **GSI** tympanometers from GSI 20 – 39 + TymStar, Allegro.



Mushroom ear tips

Part no.	Description	Unit
8518409	AKR Assortment Box, various sizes	
8500570	3-6 mm (yellow)	bag/100 pcs
8107510	4 mm (red)	bag/100 pcs
8107511	5 mm (blue)	bag/100 pcs
8500560	6 mm (green)	bag/100 pcs
8500550	7 mm (yellow)	bag/100 pcs
8500540	8 mm (red)	bag/100 pcs
8500530	9 mm (blue)	bag/100 pcs
8500520	10 mm (green)	bag/100 pcs
8500510	11 mm (yellow)	bag/100 pcs
8500500	12 mm (red)	bag/100 pcs
8500490	13 mm (blue)	bag/100 pcs
8107520	14 mm (green)	bag/100 pcs
8500480	16 mm (yellow)	bag/100 pcs
8505657	Mushroom ear tips, 5mm (blue)	bag/25 pcs
8505658	Mushroom ear tips, 6mm (green)	bag/25 pcs
8505659	Mushroom ear tips, 7mm (yellow)	bag/25 pcs
8505660	Mushroom ear tips, 8mm (red)	bag/25 pcs
8505661	Mushroom ear tips, 9mm (blue)	bag/25 pcs
8505662	Mushroom ear tips, 10mm (green)	bag/25 pcs
8505663	Mushroom ear tips, 11mm (yellow)	bag/25 pcs
8505664	Mushroom ear tips, 12mm (red)	bag/25 pcs
8505665	Mushroom ear tips, 13mm (blue)	bag/25 pcs
8505666	Mushroom ear tips, 14mm (green)	bag/25 pcs
8505667	Mushroom ear tips, 16mm (yellow)	bag/25 pcs



ABR

OAE

Save time pr. patient, earn more money and help more patients!

Use a fresh ear tip pr. patient and be sure to reach a seal right away, when inserting the ear tip. A fresh ear tip contains the optimal conditions for reaching a seal when inserting it into the patients' ear, whereas a multi-used one often will

cause you trouble and delay you in your process. By allowing yourself to plan less time pr. patient for testing, you can see more patients pr. day. That is worth much more than cleaning an ear tip, both for your clinic and your patients.

DISPOSABLE SANIBEL™ AOT and ANA ear tips



ANA ear tips for Echo-Screen®

For use with **Natus:** Echo-Screen®. **Madsen/Otometrics:** AccuScreen® Classic. **Interacoustics:** Eclipse with pink probe tip

Part no.	Description	Unit
8106363	4.0 mm (blue)	Bag/50 pcs
8106362	4.5 mm (red)	Bag/50 pcs
8106364	5.0 mm (green)	Bag/50 pcs
8106366	Infant tree (yellow)	Bag/50 pcs
8106365	Mushroom (green)	Bag/50 pcs

VOLUME DISCOUNT AVAILABLE for minimum orders of 10 bags



AOT ear tips

For use with **Madsen/Otometrics** AccuScreen®, Biologic Bio-logic® AuDX®, **Pathmedical:** all Sentiero products

Part no.	Description	Unit
8106367	3.7 mm (yellow)	Bag/50 pcs
8106368	4.0 mm (light green)	Bag/50 pcs
8106369	4.5 mm (light pink)	Bag/50 pcs
8106370	5.0 mm (light blue)	Bag/50 pcs
8106371	4-7 mm (transparent)	Bag/50 pcs

VOLUME DISCOUNT AVAILABLE for minimum orders of 10 bags



Grason & Associates ot & MO select series ear tips



OT select ear tips w/ triangle holes

Ear tips with triangle holes for use with **Otodynamics** UGD, SND, SGD probes

Part no.	Description	Unit
8102644	4.8 mm (blue)	Bag/100 pcs
8102645	6 mm (green)	Bag/100 pcs
8102646	8 mm (red)	Bag/100 pcs
8102647	10 mm (blue)	Bag/100 pcs
8102648	12 mm (green)	Bag/100 pcs
8102649	14 mm (red)	Bag/100 pcs
8102650	16 mm (blue)	Bag/100 pcs



OT select ear tips w/two probe holes

Ear tips with two probe holes for use with **Otodynamics** ECP, UGS, SNS, SGS probes.

Part no.	Description	Unit
8102632	3.5 mm (blue)	Bag/100 pcs
8102633	4.5 mm (blue)	Bag/100 pcs
8102634	5.5 mm (blue)	Bag/100 pcs
8102635	6.5 mm (blue)	Bag/100 pcs
8102636	4.8 mm (blue)	Bag/100 pcs
8102637	6 mm (green)	Bag/100 pcs
8102638	8 mm (red)	Bag/100 pcs
8102639	10 mm (blue)	Bag/100 pcs
8102640	12 mm (green)	Bag/100 pcs
8102641	14 mm (red)	Bag/100 pcs
8102642	16 mm (blue)	Bag/100 pcs



MO select ear tips for Madsen OTOflex

For **MADSEN®** OTOflex™.

Part no.	Description	Unit
8102653	3-6 mm (yellow)	Bag/100 pcs
8102654	6 mm (red)	Bag/100 pcs
8102655	7 mm (blue)	Bag/100 pcs
8102656	8 mm (green)	Bag/100 pcs
8102657	9 mm (yellow)	Bag/100 pcs
8102658	10 mm (red)	Bag/100 pcs
8102659	11 mm (blue)	Bag/100 pcs
8102660	12 mm (green)	Bag/100 pcs
8102661	13 mm (yellow)	Bag/100 pcs
8102662	14 mm (red)	Bag/100 pcs
8102663	16 mm (blue)	Bag/100 pcs
8102664	18 mm (green)	Bag/100 pcs



Ear tips for OAE screeners

For use with **Biologic AuDX**, **Welch Allyn** OAE screener and **Natus** Echo-Screen® III.



OAE ear tips, AuDX

Part no.	Description	Unit
8107885	AuDX Premie 3 mm	Bag/100 pcs
8515505	AuDX Premie 3 mm	Bag/20 pcs
8107889	AuDX Tree 3.5-5mm	Bag/100 pcs
8515506	AuDX Tree 3.5-5mm	Bag/20 pcs
8107887	AuDX Infant 6 mm	Bag/100 pcs
8515507	AuDX Infant 6 mm	Bag/20 pcs
8102681	Adult OAE foam	Bag/100 pcs
8515436	Adult OAE foam	Bag/20 pcs
8102682	Pediatric OAE foam	Bag/100 pcs
8515435	Pediatric OAE foam	Bag/20 pcs
8103123	Jumbo OAE foam	Bag/100 pcs
8515438	Jumbo OAE foam	Bag/20 pcs





We know that, first and foremost, great quality is at the top of your list, and you'll find it in each and every purchase you make from Sanibel Supply.

Bringing you the products you need is what drives us every day.

If we can ever be of assistance, don't hesitate to let us know.