

Software

Paperless ECG

seca Software enables not only the ECG trace, but also a choice of format and measurements etc. to be attached into the individual patient record. Once in the computer system the trace can be printed, faxed or emailed as desired.

- Transmission of the ECG trace and patient data through one touch function key on the ECG machine.
- User selection of viewing choices (e.g. one page report or 3 x 4, etc.).
- All lead measurement figures included in patient record.

seca Passport Patient Record Package

This software receives the trace from the ECG machine and enables attachment to the individual patient record for storage and viewing throughout the practice network.

Compatible with emis, SystmOne, vision, isoft and others.

For use with the seca CT8000i / CT8000P / CT8000L.



CT8000i Interpretive



CT8000P Interpretive



CT8000L Interpretive

seca EMR Passport Software

This software receives the trace from the ECG machine and enables attachment to the individual patient record for storage and viewing throughout the practice network.

Compatible with emis, SystmOne, vision, isoft and others.

For use with the seca CTCardioPad / CTCardioPad Mini.



seca CTCardioPad

seca CTCardioPad mini



seca CT500 ECG Trolley

- Dedicated for use with the seca range of ECG machines.
- Spring loaded arm for control of patient cable and storage basket for electrodes and paper.
- Non-metallic top.
- **3 year warranty**

