



Your cardio-respiratory partner

Data you can rely on.  
People you can trust.



## Vitalograph Peak Flow Meters

### Peak flow monitoring for adults and children

Vitalograph peak flow meters combine high quality and attractive styling with great value. These robust and easy to clean peak flow meters conform to the ISO standard for Peak Expiratory Flow Meters.

### myPEF – the only web enabled peak flow meter

myPEF™ is supplied with a unique ID card allowing users free access to a secure website where they can enter their PEF readings, eliminating the need for paper records and supporting data accuracy and effective asthma management. Suitable for use with any web enabled device.

- Easy to read graphical data display featuring data trending, personal best and asthma management plan
- Report generation featuring plotting, trending and time stamping
- Clinicians can be granted access to view patient's PEF results online
- Precise 100mm scale - longer than any other PFM's scale
- Integral dust cover for hygiene and protection
- Lightweight, pocket size design for convenient portability

### asmaPLAN

asmaPLAN™ is the first and original home use personal monitor for asthma management. Developed in the early 1970's it has evolved to become the gold standard of simple mechanical peak flow meters. The lightweight, proven and truly robust asmaPLAN is highly accurate and easy to use by both adults and children.

- |                                    |                   |
|------------------------------------|-------------------|
| ➤ Proven                           | ➤ Easy to use     |
| ➤ Low cost                         | ➤ Highly accurate |
| ➤ Suitable for adults and children | ➤ Robust          |



## Technical specification:

**Product:** Peak Flow Meter

**Model Number:** 4300

**Parameters Measured:** PEF (L/min)

**Device Dimensions:**

138mm x 56mm x 23mm (myPEF)

152mm x 61mm x 36mm (asmaPLAN)

**Nett Weight:** 79g (myPEF); 96g (asmaPLAN)

**Environmental Data:** Storage humidity: 10% to 95%, Storage temperature: 0-50°C

**Recommended Operating Temperature**

**Range:** 17 - 37°C; Design limits 10-40°C

**Ambient Humidity:** 0% to 99%

**Medical Safety Standard:** Medical Devices Directive 2007/47/EC

**Designed and Manufactured Under:** ISO

13485:2003, FDA 21 CFR 820, CMDR SOR

98/282 and JPAL

## Ordering info:

**43950** 4300 Peak Flow Meter myPEF Web Enabled

**43602** 4300 Peak Flow Meter asmaPLAN

**43703** 4300 Peak Flow Meter asmaPLAN Child

**43982** 4300 Peak Flow Meter myPEF Disposable (100)

**20242** 2024 SafeTway® Mouthpieces (200)

**29805** Mouthpiece Adaptors myPEF (96)

Vitalograph is a world leading provider of outstanding quality cardio-respiratory diagnostic devices, clinical trial services and medical equipment servicing. With a pioneering heritage of excellence spanning half a century Vitalograph continues to make valuable contributions to effective medical care and enhanced quality of life.



Vitalograph Ltd. UK  
Maids Moreton, Buckingham  
MK18 1SW  
England  
Phone: (01280) 827110  
Fax: (01280) 823302  
e-mail: sales@vitalograph.co.uk  
www.vitalograph.co.uk

Vitalograph Inc.  
13310 West 99th Street  
Lenexa, Kansas, 66215  
USA  
Toll Free: 800 255 6626  
Phone: (913) 730 3200  
Fax: (913) 730 3232  
e-mail: vitcs@vitalograph.com  
www.vitalograph.com

Vitalograph Ltd. Export  
Maids Moreton, Buckingham  
MK18 1SW  
England  
Phone: +44 1280 827120  
Fax: +44 1280 823302  
e-mail: sales@vitalograph.co.uk  
www.vitalograph.eu

Vitalograph (Ireland) Ltd.  
Gort Road Business Park  
Ennis, Co. Clare  
Ireland  
Phone: +353 65 6864100  
Fax: +353 65 6829289  
e-mail: sales@vitalograph.ie  
www.vitalograph.ie

Vitalograph GmbH  
Rellinger Straße 64a  
D-20257 Hamburg  
Germany  
Phone: (040) 547391-0  
Fax: (040) 547391-40  
e-mail: info@vitalograph.de  
www.vitalograph.de

Vitalograph Ltd. Hong Kong  
Unit 230, Corporation Park  
11 On Lai Street, Shatin  
Hong Kong  
Phone: +852 2117 2678  
Fax: +852 2117 2679  
e-mail: sales@vitalograph.hk  
www.vitalograph.hk