

3rd Edition City University Colour Vision Test and Recording Charts



38 and 24 Plate Ishihara Colour Vision Test

Colour Vision Tests

3rd Edition City University Colour Vision Test

This test displays normal protan, deutan and tritan colours, differentiation is stronger in the red/green.

Significant tritan defects may be detected if used in conjunction with the Ishihara. The test includes a screening option.

Recording Charts for City University Test

38 and 24 Plate Ishihara Colour Vision Test

Ideal for commercial/industrial colour screening. Tests for protan and deutan defects.

Patients must identify numbers or patterns formed in a matrix of small coloured dots varying in colour and intensity (pseudo-isochromatic plates).

Two editions available: 38 plate for detailed investigation – 24 plate for quick screening.

Part numbers

3rd Edition C	ity University Colour Vision Test
2206-L-1057	City University Colour Vision Test
2206-L-7125	Recording Charts for City University Colour Vision Test
38 and 24 Pla	te Ishihara Colour Vision Test
2206-P-1008	38 Plate Ishihara Colour Test
2206-P-1016	24 Plate Ishihara Colour Test